

TRULOCK™

LIGHT ACTIVATED

BAND ADHESIVE

RMO® is proud to announce a new addition of Light Activated Band Adhesive to the Trulock™ adhesive line. This new and improved formula will be replacing the existing RMO Band Cement (P/N J03280) due to recent supply issues.

The single-paste, resin-based band adhesive requires no mixing and has no working time constraints. A distinctive blue color allows you easily see excess adhesive during banding, making cleanup easy. The band adhesive will be available in either single 5g preloaded syringes or a 4 pack of 5g preloaded syringes. To place an order, please contact your RMO® Sales Representative or Customer Service.



TRULOCK® LC

TruLock LC Band Adhesive, 5g Syringe
TruLock LC Band Adhesive, 4 x 5g Syringes

PART

J04003
J04005

Ortholine B.V.

Stammerdijk 1A
1112 AA Diemen
Tel: +31 20 4165850
Fax: +31 20 4165851
Email: info@ortholine.nl
Website: www.ortholine.nl

